

Funnels

Funnels may be used for separating solids from liquids, liquids from liquids and occasionally, for pouring something into a container.

BOROSIL® brand Funnels are available in a comprehensive range for a variety of such applications. Sintered Glass Funnels are also available.


Separating Funnels come in globe, pear and cylindrical shapes with ground glass stopcocks and interchangeable stoppers. These Funnels are made of specially moulded blanks. Stopcock bores, body openings and stems are carefully aligned to facilitate the smooth and even flow of liquid with fast cut-off.



6140 - FUNNELS

Plain, 60° Angle, Short Stem


- Complies with IS 1591, DIN / ISO 4798
- Made 3.3 borosilicate glass, USP Type I
- Designed with a 60° bowl angle to ensure precise fitting of the filter paper
- High thermal shock resistance and chemical resistance
- Generally used for filtering and decanting liquids of varied densities



Product Code	Diameter mm	Quantity Per Case	Price / Piece ₹
6140057	25	30	80.00
6140058	35	30	105.00
6140065	50	20	115.00
6140069	65	20	130.00
6140071	75	20	130.00
6140077	100	40	175.00

6150 - TUBES FUNNEL TEST TUBES


- Specially Designed for Soil testing laboratory




Product Code	Capacity ml	Quantity Per Case	Price / Piece ₹
6150007	15	10	155.00
6150010	30	10	195.00

6220 - FUNNELS, POWDER

- Ideal for transferring powders
- Designed with a short wide stem to facilitate smooth flow
- Made of chemically resistant 3.3 borosilicate glass



Product Code	Capacity ml	Quantity Per Case	Price / Piece ₹
6220080	125	30	365.00
6220081	150	20	390.00
6220087	200	4	1050.00
6220089	250	4	1515.00




6230 - FUNNELS

Powder, Stem, with Cone

- Ideal for transferring powders
- Designed with a short wide stem to facilitate smooth flow
- Made of chemically resistant 3.3 borosilicate glass
- Supplied with ground glass cone for easy fitting on flask
- Joint size 24 / 29 for all capacities

Product Code	Capacity ml	Quantity Per Case	Price / Piece ₹
6230080	125	20	425.00
6230081	150	20	455.00
6230087	200	4	1215.00
6230089	250	4	1750.00



6340-FUNNELS

Separating, Globe Shape, with Stopcock and Inter- changeable Stopper

- Made from 3.3 borosilicate glass, USP Type I
- Equipped with a standard ground glass stopcock
- Globe shape, with a stem
- Generally used in the separation phase

Product Code	Capacity ml	Stopper Size	Quantity Per Case	Price / Piece ₹
6340017	125	19 / 26	10	1065.00
6340021	250	24 / 29	10	1165.00
6340024	500	24 / 29	10	1465.00
6340029	1000	29 / 32	5	1850.00
6340030	2000	29 / 32	10	3125.00
6340033	5000	34 / 35	1	7095.00



6400-FUNNELS

Separating, Pear Shape, with glass Stopcock and interchangeable Stopper

- Complies with IS 1575, DIN / ISO 4800
- Made 3.3 borosilicate glass, USP Type I
- Equipped with a standard glass stopcock
- Pear shape, with a stem
- Generally used in the separation phase

Product Code	Capacity ml	Stopper Size	Quantity Per Case	Price / Piece ₹
6400013	60*	14 / 23	10	895.00
6400017	125*	19 / 26	10	975.00
6400021	250	24 / 29	10	1080.00
6400024	500	24 / 29	10	1425.00
6400029	1000	29 / 32	10	1860.00
6400030	2000	29 / 32	10	3385.00

* Not covered in standards.

6402-FUNNELS

Separating, Pear Shape, Fitted with Boroflo Stopcock with PTFE Key and Interchangeable Stopper

- Complies with IS 1575, DIN / ISO 4800
- Made 3.3 borosilicate glass, USP Type I
- Equipped with a Boroflo stopcock with PTFE key
- Pear shape, with a stem
- Generally used in the separation phase



Product Code	Capacity ml	Stopper Size	Quantity Per Case	Price / Piece ₹
6402013	60*	14 / 23	10	800.00
6402017	125*	19 / 26	10	910.00
6402021	250	24 / 29	10	1010.00
6402024	500	24 / 29	10	1360.00
6402029	1000	29 / 32	10	1700.00
6402030	2000	29 / 32	10	2990.00

* Not covered in standards.

6403-FUNNELS

Separating, Pear Shape, with PTFE Key Stopcock, Interchangeable Stopper

- Complies with IS 1575, DIN / ISO 4800
- Made of 3.3 borosilicate glass
- USP Type I glass
- Equipped with a PTFE key stopcock
- Pear shape, with a stem
- Generally used in the separation phase



Product Code	Capacity ml	Stopper Size	Quantity Per Case	Price / Piece ₹
6403017	125*	19 / 26	10	1135.00
6403021	250	24 / 29	10	1240.00
6403024	500	24 / 29	10	1610.00
6403029	1000	29 / 32	10	1920.00
6403030	2000	29 / 32	10	3580.00

* Not covered in standards.


6405 - FUNNEL

Pressure Equilising Cylindrical, with PTFE Stopcock

- Made 3.3 borosilicate glass, USP Type I
- Equipped with totally grease free and self lubricating PTFE stopcock
- Supplied with polypropylene stopper



Product Code	Capacity ml	Socket	Cone	Quantity Per Case	Price / Piece ₹
6405012	50	19 / 26	19 / 26	5	1075.00
6405016	100	19 / 26	19 / 26	5	1285.00
6405021	250	24 / 29	24 / 29	5	1475.00
6405024	500	24 / 29	24 / 29	5	2285.00



6408 - DROPPING FUNNEL WITH GLASS STOPCOCK

- Complies with IS 1575, DIN / ISO 4800

Product Code	Capacity ml	Neck Joint	Quantity Per Case	Price / Piece ₹
6408012	50	14 / 23	5	660.00
6408016	100	19 / 26	5	720.00
6408021	250	24 / 29	5	870.00
6408024	500	24 / 29	5	1165.00
6408029	1000	29 / 32	5	1930.00




6560 - CONE Plain Shank, Single

- Made from 3.3 borosilicate glass, USP Type I
- Complies to DIN 12249
- These joints are made with precision grinding

Product Code	Cone Joint Size	Total Length	Quantity Per Case	Price / Piece ₹
6560010	10 / 19	120	5	125.00
6560012	12 / 21	120	5	135.00
6560014	14 / 23	120	5	150.00
6560019	19 / 26	125	5	185.00
6560024	24 / 29	135	5	220.00
6560029	29 / 32	135	5	275.00
6560034	34 / 35	150	5	390.00
6560040	40 / 38*	150	5	445.00
6560045	45 / 40*	150	5	530.00
6560050	50 / 42*	150	5	615.00
6560055	55 / 44	150	5	725.00

* Not covered in standard



6561- CONE Plain Shank, Double

- Made from 3.3 borosilicate glass, USP Type I
- Complies to DIN 12249
- These joints are made with precision grinding

Product Code	Cone Joint Size	Total Length	Quantity Per Case	Price / Piece ₹
6561010	10 / 19	125	5	250.00
6561012	12 / 21	135	5	280.00
6561014	14 / 23	135	5	345.00
6561019	19 / 26	140	5	390.00
6561024	24 / 29	170	5	470.00
6561029	29 / 32	170	5	640.00
6561034	34 / 35	180	5	750.00
6561040	40 / 38	185	5	640.00
6561045	45/40	185	5	750.00
6561050	50/42	195	5	800.00
6561055	55/44	195	5	925.00

6580 - SOCKETS

Plain Shank, Single

- Made from 3.3 borosilicate glass, USP Type I
- Complies to DIN 12249
- These joints are made with precision grinding



Product Code	Cone Joint Size	Total Length	Quantity Per Case	Price / Piece ₹
6580010	10 / 19	120	5	135.00
6580012	12 / 21	120	5	160.00
6580014	14 / 23	120	5	180.00
6580019	19 / 26	125	5	185.00
6580024	24 / 29	135	5	225.00
6580029	29 / 32	135	5	275.00
6580034	34 / 35	150	5	390.00
6580040	40 / 38	150	5	475.00
6580045	45 / 40	150	5	615.00
6580050	50 / 42	150	5	695.00
6580055	55 / 44	150	5	865.00

6581 - SOCKETS

Plain Shank, Double

- Made from 3.3 borosilicate glass, USP Type I
- Complies to DIN 12249
- These joints are made with precision grinding



Product Code	Cone Joint Size	Total Length	Quantity Per Case	Price / Piece ₹
6581010	10 / 19	125	5	225.00
6581012	12 / 21	135	5	270.00
6581014	14 / 23	135	5	310.00
6581019	19 / 26	140	5	350.00
6581024	24 / 29	170	5	390.00
6581029	29 / 32	170	5	470.00
6581034	34 / 35	180	5	640.00
6581040	40 / 38	185	5	695.00
6581045	45 / 40	185	5	725.00
6581050	50 / 42	195	5	775.00
6581055	55 / 44	195	5	925.00