

Cones


Cones are used to determine small amounts of suspended matter in sewage by the Imhoff Sedimentation Method.



2160 - CONES

Imhoff, Sediment, Sharp Tip

- Complies with IS 7232
- Graduated from 0 to 1 ml in 0.1 ml divisions
- Graduated from 1 to 10 ml in 0.5 ml divisions
- 10 to 40 ml in 1 ml divisions also marked at 1000 ml
- Facilitates the determination of the amount of sediment in sewage or industrial waste




Product Code	Approx Total Height mm	Approx Top Diameter mm	Quantity Per Case	Price / Piece ₹
2160094	454	110	1	4990.00

2180 - CONES

Imhoff, Sediment, Blunt Tip

- Complies with IS 7232
- The blunt tip increases its ruggedness and cleaning ease
- Ideal for use with heavy sediments where the ability to read small volumes is unimportant
- Facilitates the determination of the amount of sediment in sewage or industrial waste



Product Code	Approx Total Height	Approx Top Diameter	Quantity Per Case	Price / Piece ₹
2180093	439	110	1	5092.00