

Desiccators

BOROSIL® brand Desiccators are equipped with bodies and covers that have great strength and utility. The body is designed with a streamlined contour that facilitates easy cleaning of the used desiccant. Any desiccant such as sulphuric acid, phosphoric pentoxide, calcium chloride or silica gel may be used in Borosil brand Desiccators. For your convenience, all Borosil brand Plain Desiccators of the same size have interchangeable covers and bodies. When the ground surface is cleaned and greased, Vacuum type Desiccators hold a vacuum of 500 mm of mercury.




3082 - DESICCATORS

with Cover and Porcelain Plate, Plastic Knob

- Complies with IS 6128
- Robust Desiccator base and cover made of borosilicate glass
- Simple design for drying under atmospheric pressure
- Ground glass flanges
- Plates are positioned on an internal ledge within the base
- Generally used for drying of moist products, and for storage of moisture-sensitive substances

Product Code	Approx ID Ground Flange mm	Quantity Per Case	Price / Piece ₹
3082041	100	1	6430.00
3082042	150	1	7505.00
3082043	200	1	10465.00
3082044	250	1	13365.00
3082045	300	1	18445.00



3083 - VACUUM DESICCATORS

Vacuum, Stopcock with PTFE Spindle, and Porcelain Plate

- Complies with IS 6128
- Designed for use under vacuum
- Sturdy construction of borosilicate glass
- Ground glass flanges
- Plates are positioned on an internal ledge within the base
- Comes with glass stopcock at the top
- Generally used for drying of moist products, and for storage of moisture sensitive substances

Product Code	Approx ID Ground Flange mm	Quantity Per Case	Price / Piece ₹
3083041	100	1	9950.00
3083042	150	1	14960.00
3083043	200	1	16795.00
3083044	250	1	23315.00
3083045	300	1	29295.00

3084 - DESICCATORS

with Cover and Porcelain Plate, Plastic Knob,
Amber



- Complies with IS 6128
- Robust Desiccator base and cover made of borosilicate glass
- Simple design for drying under atmospheric pressure
- Ground glass flanges
- Plates are positioned on an internal ledge within the base
- Generally used for drying of moist products, and for storage of moisture-sensitive & light sensitive substances



Desiccators

Product Code	Approx ID Ground Flange mm	Quantity Per Case	Price / Piece ₹
3084041	100	1	8500.00
3084042	150	1	9740.00
3084043	200	1	12555.00
3084044	250	1	15950.00
3084045	300	1	24965.00

3085 - VACUUM DESICCATORS

Vacuum, stopcock with PTFE spindle and
Porcelain plate, Amber



- Complies with IS 6128
- Designed for use under vacuum
- Sturdy construction of borosilicate glass
- Ground glass flanges
- Plates are positioned on an internal ledge within the base
- Comes with glass stopcock at the top
- Generally used for drying of moist products, and for storage of moisture sensitive & light sensitive substances



Product Code	Approx ID Ground Flange mm	Quantity Per Case	Price / Piece ₹
3085041	100	1	15445.00
3085042	150	1	19840.00
3085043	200	1	21685.00
3085044	250	1	31250.00
3085045	300	1	35055.00