

# Dishes

**BOROSIL®** brand Dishes are designed to fulfil your wide-ranging needs - Crystallizing, Evaporating, Drying, Storing, Petri dishes for culture, assay and micro-biological works. All **BOROSIL®** brand Dishes are safe, durable and economical.


**BOROSIL®** brand Dishes are inert, can be sterilised and remain clear even after repeated use and washing. With each subsequent use, the cost per Dish decreases.



## 3140 - DISHES

### Crystallizing

- Reinforced and fire-polished rims reduce chipping
- Large labelling field for easy marking
- Very high chemical resistance
- High resistance to thermal shocks




Product Code	Approx O.D. x Height mm	Approx Capacity ml	Quantity Per Case	Price / Piece ₹
3140077	100 x 50	330	10	410.00
3140081	150 x 75	1150	10	715.00
3140086	190 x 100	2500	4	1185.00



## 3160 - DISHES

### Culture, Petri

- Complies with IS 2626
- These Dishes are clear and withstand repeated sterilisation (wet or dry)
- Made of 3.3 borosilicate glass using a special manufacturing process
- Beaded edges impart greater mechanical strength
- Ideal for use in routine work
- Specially formed to ensure even wall thickness and consistent optical performance



Product Code	Approx O.D. x Height mm	Quantity Per Case	Price / Piece ₹
3160065	50 x 17	100	170.00
3160072	80 x 17	100	215.00
3160077	100 x 17 <sup>#</sup>	100	255.00
3160081	150 x 20	40	625.00
3160087	200 x 20	10	1350.00


# Not covered in standard



## 3165 - DISHES

### Culture Petri, Borosil S-Line


- The inner surface of these Dishes have a very good flatness
- The accuracy of many biological assays depends on the flatness of the inside base of the Petri dish
- Made of soda lime glass ensuring that the inner surface of the top and base are absolutely flat
- These Dishes have high clarity and are free of bubbles
- Ideal for use in routine work



Product Code	Approx O.D. x Height mm	Quantity Per Case	Price / Piece ₹
3165065	50 x 12	100	70.00
3165072	80 x 15	100	80.00
3165A75	90 x 15	100	80.00
3165077	100 x 15	100	85.00
3165A77	100 x 20	100	155.00
3165081	150 x 25	40	190.00
3165087	200 x 30	10	405.00



• Dishes  
• Distilling Apparatus  
• Water Distilling Unit  
• Apparatus




## 3170 - TRAYS

### Drying, Heavy Wall

- The trays withstand hot air or steam sterilization
- Facilitates paper chromatography
- Dimensions are approximate, representing the inside length x width x height in mm

Product Code	Approx Dimension mm	Approx Capacity ml	Quantity Per Case	Price / Piece ₹
3170044	320 x 180 x 51	1500	5	<b>535.00</b>
<b>new</b> 3170A44	258 x 165 x 43	1500	5	<b>535.00</b>
3170043	350 x 200 x 58	2400	5	<b>745.00</b>
<b>new</b> 3170A43	304 x 190 x 44	2000	5	<b>745.00</b>
3170042	386 x 217 x 63	3000	5	<b>1040.00</b>
<b>new</b> 3170A42	328 x 203 x 56	3000	5	<b>1040.00</b>
3170041	404 x 257 x 61	3500	5	<b>1245.00</b>



## 3180 - Dishes

### Evaporating, Flat Bottom, With Pour Out

- Made of 3.3 borosilicate glass
- Very high chemical resistance
- High resistance to thermal shocks
- Designed with a spout to facilitate clean pouring
- Withstand to repeated sterilisation (Wet or dry)
- Beaded edge give greater mechanical strength
- Ideal for use in routine work


Product Code	Approx O.D. x Height mm	Approx Capacity ml	Quantity Per Case	Price / Piece ₹
3180072	80 x 45	165	10	<b>350.00</b>
3180078	105 x 55	290	10	<b>495.00</b>
3180081	150 x 80	855	10	<b>925.00</b>
3180087	200 x 100	1790	8	<b>1705.00</b>

# BOROSIL®



On line  
Certificate  
available on  
our Website  
[www.borosil.com](http://www.borosil.com)


Individually Calibrated  
**Volumetric Flasks**  
with actual meniscus photograph



## State of the art calibration

- Laser identification of meniscus
- Actual photo of individual meniscus captured on certificates (for 5ml to 100 ml)

Now ACCESS your  
BOROSIL glassware certification  
from anywhere, anytime.....



Download the **BOROSIL®** app

