

BOROSIL® Lab Accessories

TONGS FOR BEAKERS & FLASKS

- Made of stainless steel
- Used to handle hot beakers & flasks
- Provided with silicon sleeves or in fibre glass yarn to resist direct flame



Product Code	Description	Size Inches	Case Lot No.	Price / Piece ₹
LATG8888P10	Beaker tong	10	6	324.00
LATG8888B10	Beaker tong	12	6	427.00
LATG8888F10	Flask tong	10	6	407.00
new LATG8888F12	Flask tong	12	6	510.00



TONGS FOR CRUCIBLE

- Made of Stainless Steel
- Used to handle hot crucibles

Product Code	Description	Size inches	Case Lot No.	Price / Piece ₹
LACT8888006	Crucible tong with bow stainless steel	6	6	165.00
LACT8888008	Crucible tong with bow stainless steel	8	6	206.00
new LACT8888010	Tong crucible, SS,	10	6	280.00
LACT8888012	Crucible tong with bow stainless steel	12	6	376.00



CLAMPS

- Available in Zinc & Aluminium diecasted
- It is coated with LEAD-FREE powder coating & fitted with brass screw which makes it totally RUST-FREE CLAMP
- Two & three finger clamp Jaws can be covered with PVC sleeves
- Two & four finger clamp has cork lined Jaws holds upto 90 mm articles
- It is provided with 8 mm diameter & 150 mm long Stainless Steel(Grade 201)rod

Product	Description	Case Lot No.	Price / Piece ₹
LAFC8888Z02	Two finger clamp, Zinc pressure die casted	10	324.00
LAFC8888AD3	Three finger clamp, Aluminum pressure die casted	10	376.00
LAFC8888ZD4	Four finger clamp, Zinc pressure die casted	10	324.00
LAFC8888AD4	Four finger clamp Aluminum pressure die casted	10	345.00

CLAMPS WITH BOSSHEAD

- Available in Zinc & Aluminium diecasted
- It is coated with LEAD-FREE powder coating & fitted with brass screw which makes it totally RUST-FREE CLAMP
- Two & three finger clamp Jaws can be covered with PVC sleeves
- Two & four finger clamp has cork lined Jaws holds upto 90 mm articles
- With adjustable bosshead attached clamp length can be vary, can be rotated 360° and locked at any angle

Product Code	Description	Case Lot No.	Price / Piece ₹
L AFC8888Z03	Two finger clamp (Zinc PDC) with adjustable bosshead	12	587.00
L AFC8888A03	Three finger clamp (Aluminium PDC) with adjustable bosshead	12	587.00
L AFC8888Z04	Four finger clamp (Zinc PDC) with adjustable bosshead	12	644.00
L AFC8888A04	Four finger clamp (Aluminium PDC) with adjustable bosshead	12	670.00



UNIVERSAL CLAMPS

- Single Burette / Test Tube Clamp made of stainless steel Jaws & Aluminium bosshead
- Round jaws are covered with PVC coating
- Opens from 15 mm to 45 mm in Diameter
- Clamps attaches to rods up to 13 mm in diameter. Threaded rod and wing nut adjust clamp angle and overall length from 13 cm to 15 cm


Product Code	Description	Case Lot No.	Price / Piece ₹
L ABC8888SC1	Single Burette clamp, PVC coated jaws stainless steel	12	427.00
L ABC8888S01	Single Burette clamp, PVC uncoated jaws, stainless steel	12	376.00




BURETTE CLAMPS

- Single / Double Burette clamp is made of Aluminium die casting
- Feature black PVC sleeves for firm gripping of burettes.
- This adjustable burette clamp is suitable for all burettes up to 100 ml which are held in place by the pressure created by the springs yet allows easy operation without damaging glassware.

Product Code	Description	Case Lot No.	Price / Piece ₹
L ABC8888B01	Single burette clamp, Aluminium PDC Fig. A	6	644.00
L ABC8888A02	Double burette clamp, Aluminium PDC Fig. B	6	644.00
L ABC8888AE1	Single burette clamp, Aluminium PDC (Eco.) Fig. C	6	324.00
L ABC8888AE2	Double burette clamp, Aluminium PDC (Eco.) Fig. D	6	427.00





BOSSHEDS

- Bosshead made of Zinc & Aluminum pressure die casted
- It can hold a rod upto 16 mm diameter at 90° angle and Swivel bosshead rotate at 360° angle
- It is coated with LEAD-FREE powder coating & fitted with Nickel plated thumb screw

Product Code	Description	Case Lot No.	Price / Piece ₹
LABH8888RBL	Regular bosshead Zinc PDC Fig. A	12	165.00
LABH8888OBL	Round bosshead Zinc PDC Fig. B	12	206.00
LABH8888JBL	Jumbo bosshead Aluminum PDC Fig C	12	427.00
LABH8888SBL	Swivel bosshead Zinc PDC Fig. D	12	288.00

RETORT RING

- Strong circular mild steel ring with 6 mm diameter stem
- Open bosshead fitted with thumb screw to accept rods from 8 mm to 15 mm diameter



Product Code	Description	Size in inches	Case Lot No.	Price / Piece ₹
LARR8888003	Retort Ring with bosshead	3	6	268.00
LARR8888004	Retort Ring with bosshead	4	6	335.00

RETORT BASE

- Superior Heavy duty steel alloy stand with rubber lining on a shorter side gives extra stability
- With lead-free powder coating



Product Code	Description	Size inches	Case Lot No.	Price / Piece ₹
LARB8888H06	Retort base, Heavy	6" x 4"	6	587.00
LARB8888H08	Retort base, Heavy	8" x 5"	6	695.00
LARB8888H10	Retort base, Heavy	10" x 6"	6	1174.00

ROD FOR RETORT BASE

- With stainless steel Grade 201



Product Code	Description	Approx O.D. x Height mm	Case Lot No.	Price / Piece ₹
LARS8888500	Retort Rod (SS Grade 201)	12 x 500	2	407.00
LARS8888600	Retort Rod (SS Grade 201)	12 x 600	2	494.00
LARS8888750	Retort Rod (SS Grade 201)	12 x 750	2	567.00

TEST TUBE STAND

- It is made of Aluminum sheet, Z type, new design, single piece fabrication
- One row of holes of 30 mm diameter, the second row of holes of 25 mm diameter and third row of hole of 16 mm diameter
- This aluminium Z shape design is Coated with lead free powder coating make it corrosion resistant

Product Code	Description	Case Lot No.	Price / Piece ₹
LATS8888H19	Z' Shape Test Tube Stand, Aluminum, 19 hole assorted	6	427.00
LATS8888H25	'Z' Shape Test tube stand, 25 hole for 16mm tube	6	427.00



TEST TUBE HOLDER

- A pair of spring strip of hard rolled stainless steel, mounted in wooden handle in cross-scissor shape, and accepts tubes of diameter 12-43 mm.
- Overall length about 225mm


Product Code	Description	Case Lot No.	Price / Piece ₹
LATH8888000	Test tube holder, cross pattern, Stainless steel	24	185.00



BURNERS


- Bunsen Burner 13mm/11mm ideal for use with natural gas and L.P.G.
- Nickel-plated brass burner tube with Rotatable air regulator
- Enamelled Zinc pressure die casted base
- Base diameter 80 mm

Product Code	Description	Caselot No.	Price / Piece ₹
LABR8888Z11	Bunsen burner, Zinc PDC base with 11mm brass pipe for L.P. G. Fig. A	2	803.00
LABR8888Z13	Bunsen burner, Zinc PDC base with 13mm brass pipe for L.P. G. Fig. B	2	803.00
LABR8888013	Burner, sheet metal base with 13mm pipe, for L.P. G. Fig. C	2	644.00



SPATULAS

- Made of stainless steel grade 202
- Spatulas are used for measuring solids moving objects and scraping material out of
- beakers



Product Code	Description	Size in inches	Case Lot No.	Price / Piece ₹
LACS8888006	Chattaway Spatula	6	20	67.00
LACS8888008	Chattaway Spatula	8	20	77.00
LASL8888C08	Micro Chattaway Spatula	8	20	72.00
LASC8888M06	Micro Scoop	6	10	57.00
LASC8888SM6	Semi Micro Spatula	6	10	82.00
LASL8888075	Big Spoon Spatula	7.5	10	124.00
LASL8888009	Big Spoon Spatula	9	10	139.00
LASP8888B06	Both side flat	6	20	72.00
LASP8888B08	Both side flat	8	20	93.00
LASP8888004	One end flat and one end Spoon	4	50	52.00
LASP8888006	One end flat and one end Spoon	6	50	57.00
LASP8888008	One end flat and one end Spoon	8	50	82.00
LASP8888010	One end flat and one end Spoon	10	50	129.00
LASP8888012	One end flat and one end Spoon	12	50	149.00
LASO8888W04	Ointment Spatula	4	10	165.00
LASO8888W06	Ointment Spatula Wooden Handle	6	10	175.00
LATG8888S06	Ointment Spatula Forged Stainless Steel	6	10	350.00
LASL8888P06	Pointer Spatula	6	10	98.00
LASC8888165	165 mm Scoop	6.5	10	82.00
LASL8888510	Big Spoon Spatula	10	10	242.00

* For products listed on the above table use HSN code 90189022

FORCEPS POINTED

Product Code	Description	Case Lot No.	Price / Piece ₹
LAFP8888004	Forceps, Pointed dissecting, 4"	10	100.00
LAFP8888005	Forceps, Pointed dissecting, 5"	10	130.00
LAFP8888006	Forceps, Pointed dissecting, 6"	10	160.00

FORCEPS PTFE POINTED

Product Code	Description	Case Lot No.	Price / Piece ₹
LAFPP8888004	Forceps, <u>PTFE</u> , Pointed end, 4"	10	1400.00
LAFPP8888006	Forceps, <u>PTFE</u> , Pointed end, 6"	10	1800.00

FORCEPS BLUNT DISSECTING

Product Code	Description	Case Lot No.	Price / Piece ₹
LAFB8888005	Forceps, Blunt dissecting, 5"	10	140.00
LAFB8888006	Forceps, Blunt dissecting, 6"	10	150.00
LAFB8888008	Forceps, Blunt dissecting, 8"	10	210.00

FORCEPS PTFE BLUNT

Product Code	Description	Case Lot No.	Price / Piece ₹
LAFBP8888004	Forceps, PTFE, Blunt end, 4"	10	1400.00
LAFBP8888006	Forceps, <u>PTFE</u> , Blunt end, 6"	10	1800.00

WIRE GAUGE

Product Code	Description	Case Lot No.	Price / Piece ₹
LAWG8888C06	Wire Gauge with Ceramic center 5" x 5"	10	120.00
LAWG8888C05	Wire Gauge with Ceramic center 6" x 6"	10	95.00

TRIANGULAR TRIPOD

Product Code	Description	Case Lot No.	Price / Piece ₹
LATT8888200	Triangular Tripod, SS, 200 x 125 mm	2	560.00

new

new

new

new

new

new