

Premium Micro Pipette Tips

- Certified DNase, RNase and Pyrogen Free
- Free of PCR inhibitors, lubricants and heavy metals
- Produced from special diamond polish mould using medical grade virgin polypropylene which provides very high transparency
- With graduations for visual control of liquid transfer
- Available in both standard as well as low retention
- With no secondary treatment and with no silicone additives for low retention
- With special PE filters for prevention of aerosols and contamination
- Specifically configured to fit major Pipette Brands available worldwide

Filter tip features can be selected as per permutation of

- Sterile
- Non sterile
- Low retention
- Filter

FILTER TIPS

BOROSIL filter tips are manufactured using specifically developed Polyethylene filters to block aerosols from entering into the pipette.

This helps in prevention of cross contamination.

Ideal for use in applications like sequencing, DNA amplification, Tissue culture, RT-PCR and other critical protocols.

LOW RETENTION TIPS

BOROSIL low retention tips are manufactured from a specific resin in combination with high quality diamond finish moulds to give super hydrophobic internal surface. This results in tips having minimum retention with maximum accuracy.

Further since no silicone treatment is given, there is absolutely no leaching effect.



UNIVERSAL GRADUATED TIPS

Universal Graduated Tips Bulk Pack

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT0000010B000	BGT-10B	10		1000 tips/ pack	500.00
BGT0000010XLB0	BGT-10XLB	10	Long	1000 tips/ pack	500.00
BGT0000200B000	BGT-200B	200		1000 tips/ pack	460.00
BGT000200BC000	BGT-200BC	200	Yellow	1000 tips/ pack	460.00
BGT0001000B000	BGT-1000B	1000		1000 tips/ pack	700.00
BGT001000BC000	BGT-1000BC	1000	Blue	1000 tips/ pack	700.00
BGT001000XLB00	BGT-1000XLB	1000	Long	500 tips/ pack	840.00
BGT0C10005000B	BGT-C10005	5000	Only C1	100 tips/ pack	1140.00
BGT0C30005000B	BGT-C30005	5000	Only C3	100 tips/ pack	1140.00
BGT0C10010000B	BGT-C10010	10000	Only C1	100 tips/ pack	2500.00
BGT0C30010000B	BGT-C30010	10000	Only C3	100 tips/ pack	2500.00

Universal Graduated Tips Bulk Pack Sterile

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010SB000	BGT-10SB10	10		1000 tips/ pack	670.00
BGT000010XLSB0	BGT-10XLSB	10	Long	1000 tips/ pack	670.00
BGT000200SB000	BGT-200SB	200		1000 tips/ pack	670.00
BGT001000SB000	BGT-1000SB	1000		1000 tips/ pack	1000.00
BGT001000XLSB0	BGT-1000XLSB	1000	Long	500 tips/ pack	1230.00

Universal Graduated Tips Bulk Pack Sterile Low Retention

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010SLB000	BGT-10SLB10	10		1000 tips/ pack	920.00
BGT00010XLSLB0	BGT-10XLSLB	10	Long	1000 tips/ pack	920.00
BGT00200SLB000	BGT-200SLB200	200		1000 tips/ pack	920.00
BGT01000SLB000	BGT-1000SLB	1000		1000 tips/ pack	1500.00
BGT01000XLSLB0	BGT-1000XLSLB	1000	Long	500 tips/ pack	1840.00

Universal Graduated Tip Racks

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT0000010R000	BGT-10R10	10		10 X 96 Tips/rack	2130.00
BGT0000010XLR0	BGT-10XLR	10	Long	10 X 96 Tips/rack	2130.00
BGT0000200R000	BGT-200R	200		10 X 96 Tips/rack	2010.00
BGT0001000R000	BGT-1000R	1000		6 X 96 Tips/ rack	1290.00
BGT0001000XLR0	BGT-1000XLR	1000	Long	6 X 96 Tips/ rack	1650.00

Universal Graduated Tip Racks Sterile

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010SR000	BGT-10SR10	10		10 X 96 Tips/rack	2290.00
BGT0010XLSR000	BGT-10XLSR	10	Long	10 X 96 Tips/rack	2290.00
BGT000200SR000	BGT-200SR	200		10 X 96 Tips/rack	2180.00
BGT001000SR000	BGT-1000SR	1000		6 X 96 Tips/ rack	1510.00
BGT001000XLSR0	BGT-1000XLSR	1000	Long	6 X 96 Tips/ rack	1930.00

UNIVERSAL GRADUATED TIPS

Universal Graduated Tip Racks Sterile Low Retention

Product Code	Short Code	Tip Vol µl	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010SLR000	BGT-10SLR	10		10 X 96 Tips/rack	2450.00
BGT00010XLSLR0	BGT-10XSLR	10	Long	10 X 96 Tips/rack	2450.00
BGT00200SLR000	BGT-200SLR	200		10 X 96 Tips/rack	2400.00
BGT01000SLR000	BGT-1000SLR	1000		6 X 96 Tips/ rack	1950.00
BGT01000XLSLR0	BGT-1000XSLR	1000	Long	6 X 96 Tips/ rack	2310.00

Universal Graduated Tip Reload Stack

Product Code	Short Code	Tip Vol µl	Special Feature	Selling Unit	Price/ Unit ₹
BGT0000010Z000	BGT-10Z	10		10 X 96 tips / rack	1270.00
BGT0000200Z000	BGT-200Z	200		10 X 96 tips / rack	1270.00
BGT0001000Z000	BGT-1000Z	1000		5 X 96 tips / rack	800.00

Universal Graduated Tip Reload Stack Sterile

Product Code	Short Code	Tip Vol µl	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010SZ000	BGT-10SZ	10		10 X 96 tips / rack	1570.00
BGT000200SZ000	BGT-200SZ	200		10 X 96 tips / rack	1500.00
BGT001000SZ000	BGT-1000SZ	1000		5 x 96 tips/ rack	1120.00

Universal Graduated Tip Reload Stack Sterile Low Retention

Product Code	Short Code	Tip Vol µl	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010SLZ000	BGT-10SLZ	10		10 X 96 tips / rack	1840.00
BGT00200SLZ000	BGT-200SLZ	200		10 X 96 tips / rack	1790.00
BGT01000SLZ000	BGT-1000SLZ	1000		5 x 96 tips/ rack	1400.00

UNIVERSAL GRADUATED FILTER TIPS

Universal Graduated Filter Tips Bulk Pack

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010FB000	BGT-10FB	10		1000 tips/ pack	2290.00
BGT000010XLFB0	BGT-10XFB	10	Long	1000 tips/ pack	2290.00
BGT000020FB000	BGT-20FB	20		1000 tips/ pack	2340.00
BGT000050FB000	BGT-50FB	50		1000 tips/ pack	2340.00
BGT000100FB000	BGT-100FB	100		1000 tips/ pack	2340.00
BGT000200FB000	BGT-200FB	200		1000 tips/ pack	2340.00
BGT001000FB000	BGT-1000FB	1000		1000 tips/ pack	2900.00
BGT001000XLFB0	BGT-1000XFB	1000	Long	500 tips/ pack	3680.00

Universal Graduated Filter Tips Bulk Pack Sterile

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010SFB000	BGT-10SFB	10		1000 tips	2660.00
BGT00010XLSFB0	BGT-10XSFB	10	Long	1000 tips	2660.00
BGT00020SFB000	BGT-20SFB	20		1000 tips	2660.00
BGT00050SFB000	BGT-50SFB	50		1000 tips	2660.00
BGT00100SFB000	BGT-100SFB	100		1000 tips	2660.00
BGT00200SFB000	BGT-200SFB	200		1000 tips	2660.00
BGT01000SFB000	BGT-1000SFB	1000		1000 tips	3250.00
BGT01000XLSFB0	BGT-1000XSFB	1000	Long	500 tips	3990.00

Universal Graduated Filter Tips Bulk Pack Sterile Low Retention

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010SLFB0	BGT-10LFRB	10		1000 Tips	2830.00
BGT00010XLSLFB	BGT-10XLFR	10	Long	1000 Tips	2830.00
BGT000020SLFB0	BGT-20LFR	20		1000 Tips	2830.00
BGT000050SLFB0	BGT-50LFR	50		1000 Tips	2830.00
BGT000100SLFB0	BGT-100LFR	100		1000 Tips	2830.00
BGT000200SLFB0	BGT-200LFR	200		1000 Tips	2830.00
BGT001000SLFB0	BGT-1000LFR	1000		1000 Tips	3740.00
BGT01000XLSLFB	BGT-1000XLFR	1000	Long	500 Tips	4600.00

Universal Graduated Filter Tip Racks

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000010FR000	BGT-10FR	10		10 X 96 Tips/rack	3350.00
BGT000010XLFR0	BGT-10XFR	10	Long	10 X 96 Tips/rack	3160.00
BGT000020FR000	BGT-20FR	20		10 X 96 Tips/rack	3350.00
BGT000050FR000	BGT-50FR	50		10 X 96 Tips/rack	3350.00
BGT000100FR000	BGT-100FR	100		10 X 96 Tips/rack	3350.00
BGT000200FR000	BGT-200FR	200		10 X 96 Tips/rack	3350.00
BGT001000FR000	BGT-1000FR	1000		6 X 96 Tips/ rack	2290.00
BGT001000XLFR0	BGT-1000XFR	1000	Long	6 X 96 Tips/ rack	2910.00

UNIVERSAL GRADUATED FILTER TIPS

Universal Graduated Filter Tip Rack Sterile

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010SFR000	BGT-10SFR10	10		10 X 96 Tips/rack	3460.00
BGT00010XLSFR0	BGT-10XLSFR	10	Long	10 X 96 Tips/rack	4240.00
BGT00020SFR000	BGT-20SFR	20		10 X 96 Tips/rack	3460.00
BGT00050SFR000	BGT-50SFR	50		10 X 96 Tips/rack	3460.00
BGT00100SFR000	BGT-100SFR	100		10 X 96 Tips/rack	3460.00
BGT00200SFR000	BGT-200SFR	200		10 X 96 Tips/rack	3460.00
BGT01000SFR000	BGT-1000SFR	1000		6 X 96 Tips/ rack	2570.00
BGT01000XLSFR0	BGT-1000XSFR	1000	Long	6 X 96 Tips/ rack	3190.00

Universal Graduated Filter Tip Rack Sterile Low Retention

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT0010SLFR000	BGT-10SLFR	10		10 X 96 Tips/rack	4240.00
BGT0010XLSLFR0	BGT-10XLSLFR	10	Long	10 X 96 Tips/rack	4240.00
BGT0020SLFR000	BGT-20SLFR	20		10 X 96 Tips/rack	4240.00
BGT0050SLFR000	BGT-50SLFR	50		10 X 96 Tips/rack	4240.00
BGT0100SLFR000	BGT-100SLFR	100		10 X 96 Tips/rack	4240.00
BGT0200SLFR000	BGT-200SLFR	200		10 X 96 Tips/rack	4240.00
BGT1000SLFR000	BGT-1000SLFR	1000		6 X 96 Tips/ rack	3180.00
BGT1000XLSLFR0	BGT-1000XSLFR	1000	Long	6 X 96 Tips/ rack	3570.00

UNIVERSAL GRADUATED INDIVIDUALLY WRAPPED TIPS

Universal Graduated Tips Sterile Individually Wrapped

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT000200SW000	BGT-2000SW	200		400 Tips/ pack	2440.00
BGT001000SW000	BGT-1000SW	1000		200 tips/ pack	1630.00

Universal Graduated Tips Sterile Low Retention Individually Wrapped

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT00200SLW000	BGT-2000SLW	200		400 tips/ pack	2570.00
BGT01000SLW000	BGT-1000SLW	1000		200 tips/ pack	1790.00

EMPTY TIP BOX RACKS

Product Code	Short Code	Tip Vol μ l	Special Feature	Selling Unit	Price/ Unit ₹
BGT00010ETR000	BGT-10ETR	10		10 Rack Boxes	1650.00
BGT00010XLETRO	BGT-10XETR	10	Long	10 Rack Boxes	1650.00
BGT00200ETR000	BGT-200ETR	200		10 Rack Boxes	1650.00
BGT01000ETR000	BGT-1000ETR	1000		6 Rack Boxes	1100.00
BGT01000XLETRO	BGT-1000XETR	1000	Long	6 Rack Boxes	1140.00

Bulk pack is a pack of 1000 tips

Racks are tip boxes with 96 tips in a rack

Reload racks are racks of 96 tips

HSN CODE : 90319000



Micro Centrifuge Tubes

Certified RNase, DNase, DNA & Pyrogen Free | Centrifuge upto 20000 RCF
Temperature range: -80°C to 121 C | FDA approved resins | Universal rotor compatibility.

Product Code	Description	Pack Size	Price / Pack ₹
BGPMCT015C	1.5 ml Micro Centrifuge Tubes With Cap	500 x 2 (1000)	850.00

* For Products listed on the above table use HSN code 39269099

We, **BOROSIL LIMITED** (which was a wholly owned subsidiary of **Borosil Glass Works Limited** and formerly known as Hopewell Tableware Limited), a company registered under the Companies Act, 1956 and having its registered office at 1101, Crescenzo House, 44, Dr. R.G. Thadani Marg, Worli, Mumbai - 400018. are the owners and proprietor's of the Trademark **BOROSIL®**, which is a registered mark. Pursuant to NCLT Mumbai Bench order dated 15 Jan 2020, sanctioning a Composite Scheme of Amalgamation and Arrangement, the trademark **BOROSIL** which was earlier owned by BGWL has been vested in **BOROSIL LIMITED**.

The said, trademark **BOROSIL®** has been used in the market, to sell and export a wide variety of glassware products since more than 5(five) decades throughout India as well as several countries across the globe.

The Trademark **BOROSIL®** has acquired a very high reputation and goodwill by virtue of long, extensive and continuous use. The said trademark is associated and identified exclusively with the goods manufactured, marketed and sold by **BOROSIL LIMITED** (which was a wholly owned subsidiary of **Borosil Glass Works Limited** and formerly known as Hopewell Tableware Limited).

It has, from time to time come to our notice that some persons occasionally try to sell spurious goods by copying the **BOROSIL®** trademark and/or 'Borosil Logo' 

The trade and public are cautioned that only products manufactured, marketed and sold by erstwhile **Borosil Glass Works Limited** and **BOROSIL LIMITED** are permitted to carry the trademark **BOROSIL®** or 'Borosil Logo' 

Those who are deceiving the public by copying the Trademark **BOROSIL®** or 'Borosil Logo'  or using it wrongly on products not made, marketed and /or sold by **Borosil Limited** (which was a wholly owned subsidiary of **Borosil Glass Works Limited** and formerly known as Hopewell Tableware Limited), will be prosecuted and **Borosil Limited** will adopt all legal recourse available to it, which may be Civil or Criminal in nature and may include police raids, to prevent such unauthorized use without giving further notice.

The Trade and Public are cautioned that they should always request for a bill on purchase of their products, as spurious goods are harmful and have dangerous and adverse effects to the life of the user as well as the nation.

We request you to bring to our attention any instances of such nefarious activities, which you may come across, we shall maintain the secrecy of you identify.

Correspondence in this connection may please be made at registered office address of the Company.